



LOUPES

Magnification Systems



EOS

wireless

The total absence of cables, a lightweight body combined with a powerful light and a long lasting battery, this is EOS Wireless in a nutshell. The new illuminator is equipped with Advanced Hand Gesture Recognition technology.

EOS Wireless c'est l'absence totale de câbles combinée à la puissance, la légèreté et l'autonomie des batteries. Le nouvel éclairage intègre la technologie Advanced Hand Gesture Recognition.

Aseptic control

works on both sides



LEARN MORE



univetloupes.com/eos-gesture

Hand Gesture

RECOGNITION TECHNOLOGY

The - Advanced Hand Gesture Recognition - technology, working on both sides, allows you to control the device without touch to improve hygiene.

La technologie - Advanced Hand Gesture Recognition - permet la gestion des commandes mains libres pour des performances de plus en plus élevées.



ASEPTIC CONTROL
Contrôle aseptique



TOUCH CONTROL
Contrôle sans contact



SAVE BATTERY
Économiser la batterie



headlight

DETAILS

EOS
wireless

FEATURES

Advanced Hand Gesture Recognition technology
Blue light compliant / Anti-lumière bleue
Color temperature / Température de couleur: 5700 K
Brightness / Luminosité: 38000 lux
Battery life / Autonomie de la batterie: 140 min
Headlight weight / Poids lampe frontale: 36 g

SPOT SHAPE



Working distance (cm)	Ø (cm)	lux
35	7,5	38000

ACCESSORIES INCLUDED

2 Batteries
Battery charger
Universal magnetic adapter
UV 525 nm composite filter

OPTIONAL

Custom magnetic adapter (only for Univet TTL Loupes)

EOS
Wireless



EOS
Defined
light spot
Headlight System

www.univetloupes.com



UNIVET LOUPES Spa Via Giovanni Prati, 87 - 25086 Rezzato (BS) Italy +39 030 2499411 · info@univetloupes.com